


DONJOY®

QUICK REFERENCE GUIDE

FIND YOUR PERFECT BRACE!

BACK



WRIST/THUMB




ELBOW



KNEE



ANKLE



For more information please contact Pharmabroker Sales Pty Ltd on: (02) 8878 9763

Donjoy is registered trademark of DJO. DJO Australia PTY LTD, ABN 64 094 431 473, PO Box 6057, Frenchs Forest DC, NSW 2086

CONDITION GUIDE

Ankle	Donjoy Fortilax	Donjoy Strapping	Donjoy Stabilizing Brace	Aircast A60
General soreness and pain	✓			
Sprain & Strains	✓	✓	✓	✓
Swelling	✓	✓		
Mild tendonitis	✓	✓		✓
Mild osteo arthritis	✓	✓		
Proprioception	✓	✓	✓	✓
Resumption of activity	✓	✓	✓	✓
Post operative support	✓	✓	✓	✓
Resumption of sport activities			✓	✓
Protective and preventative use for sports activities			✓	✓
Chronic Laxity		✓	✓	✓

Knee	Donjoy Fortilax	Donjoy Rotulax open	Donjoy Stabilax
General soreness and pain	✓	✓	✓
Mild sprain	✓	✓	✓
Moderate sprain			✓
Mild meniscus injuries		✓	✓
Swelling	✓	✓	✓
Mild osteoarthritis	✓	✓	✓
Proprioception	✓	✓	✓
Chronic Laxity			✓
Resumption of activity		✓	✓
Resumption of sport activities		✓	✓
Patella stabilisation		✓	✓
Patellar subluxation		✓	✓
Patellar tendinitis		✓	✓
Patellofemoral syndrome		✓	✓
Post patellar operation		✓	✓


Back	Donjoy Conforstrap	Donjoy Porostrap	Donjoy Immostrap
Lower back pain	✓	✓	✓
Compression	✓	✓	✓
Postural reminder	✓	✓	✓
Abdominal support	✓	✓	✓
Muscle relaxant	✓		
Post operative support	✓	✓	✓

Wrist and thumb	Donjoy Quick-Fit (One Size) (Wrist Splint)	Donjoy THUMB-OPRENE™ (One size) (Wrist & Thumb)	Donjoy ErgoForm (Thumb Splint)
General soreness and pain	✓	✓	✓
Sprain	✓	✓	✓
Swelling	✓	✓	
Mild tendonitis	✓	✓	✓
Mild osteo arthritis	✓	✓	✓
Ligament instability	✓	✓	✓
Immobilisation	✓	✓	✓
Rehabilitation	✓	✓	✓
Carpal Tunnel syndrome	✓		
De Quervain		✓	✓
Hyperextension protection		✓	✓
Post operative support	✓	✓	✓

Elbow	Aircast armband (One size)	Donjoy Epiforce
General soreness and pain		✓
Epicondilitis (tennis & golf elbow)	✓	✓
Post operative support		✓
Bursitis		✓

SIZE CHART

Body Part	X-Small	Small	Medium	Large	X-Large	XX-Large	How to measure
Back	56-68	68-82	82-98	98-116	116 - 136		Measure around the waist at the top of the navel.
Knee	35-38	39-42	43-46	47-49	50-53	54-57	15 cm above mid knee cap/patella
Ankle*	19-20	21-22	23-24	25-26	27-28	29-30	Measure around the ankle just above the bone
Stabilizing Ankle	18-20	21-24	25-28	29-31	32-35		Measurement taken at circumference from the base of the heel around the ankle, while wearing an athletic sock.
A60	Men	up to 6	6.5-10.5	11+			Shoe size
	Women	up to 8.5	9-13	13.5+			
Elbow		22.2-24.8	24.8-27.9	27.9-31.1	31.1 - 34.3		Measure around elbow for correct sizing
Wrist	11-14	14-16	16-19	19-21	21-24		Measure around the wrist at the smallest point
Thumb			14-18	18-22			Measure around the wrist at the smallest point



*Sizing for Donjoy Fortulax and Strapping ankle braces